



## Spectrum Accessories

Shanghai Goptica Co.,Ltd

CHINA · SHANGHAI

[www.goptica.com](http://www.goptica.com)

## Catalogue

|                                                        |    |   |
|--------------------------------------------------------|----|---|
| I、 Standard white board.....                           | 1  | - |
| II、 Diffuse Reflection Grayscale Step Plate.....       | 2  | - |
| III、 Standard Aluminum Mirror.....                     | 3  | - |
| IV、 Deuterium Halide Tungsten Light Source.....        | 4  | - |
| V、 Halogen Light Source.....                           | 5  | - |
| VI、 Standard Light Source Kit.....                     | 6  | - |
| VII、 Deuterium Light Source.....                       | 8  | - |
| VIII、 Mercury-argon Light Source.....                  | 9  | - |
| IX、 Laser.....                                         | 9  | - |
| X、 Laser-driven plasma wide-spectrum light source..... | 10 | - |
| XI、 LED Light Source.....                              | 12 | - |
| XII、 Optical Holders AO Series.....                    | 13 | - |
| XIII、 Fiber Optic Collimator.....                      | 15 | - |

## Spectrum Accessories

### Standard White Board



#### Product Description

A whiteboard is a shortened term for a standard whiteboard or reference whiteboard, which is an optical transmission standard. It is a whiteboard that has been calibrated and provides spectral reflectance data, supplied by the National Metrology Institute (National Standards Administration). Broad users rely on this data to measure the optical properties of other samples. Regardless of the instrument or method used, all measurement data must ultimately be traced back to the standard of the National Metrology Institute. Only in this way can measurements from different sources be comparable, effective, accurate, and credible.


#### Optical Performance

- High reflectance ratio, requiring the reflectance ratio to be as high as possible, with a reflectance better than 98%.
- Good diffusing performance, surface without gloss, close to Lambertian diffuser.
- Low spectral selectivity, i.e., good neutrality, with the same high spectral reflectance ratio for all wavelengths, using wavelength 250-2500nm.
- Non-transparency to avoid background and edge loss.
- Good flatness and uniformity, with the entire surface flat and the same reflectance ratio at all points.
- Sufficiently good optical stability.
- The value is traceable to PRD, with its reflectance ratio known, providing the reflectance data report from the Chinese Academy of Metrology.

### Specifications

| Product name                | Standard whiteboard                                                                 |
|-----------------------------|-------------------------------------------------------------------------------------|
| Specification model         | DR-C50                                                                              |
| Whiteboard diameter for use | 50mm (10mm-3000mm customizable)                                                     |
| Whiteboard material         | Domestic PTFE, Imported PTFE (optional)                                             |
| Spectral range              | 205nm-2500nm                                                                        |
| Reflectivity                | Visible light range better than 98%                                                 |
| Data and Information        | Provide a copy of the reflectance data report from the Chinese Academy of Metrology |
| Shell material\color        | Aluminum alloy\black, blue                                                          |
| Shell size                  | 59*17mm                                                                             |
| Weight                      | 92g                                                                                 |

## Diffuse Reflection Grayscale Step Plate



### Product Description

A standard diffuse gray reference plate is independently developed by our company, processed through special techniques, and assembled to achieve high diffuse reflection, precise measurement, simple and convenient use, and easy portability.

Large-area dark, gray, and white target plates are ideal for accurately evaluating short-range and long-range sensitivity within the dynamic range of LiDAR systems.

### Optical Performance

- The reflective surface of the gray reference plate is an ideal Lambertian diffuse surface, with reflected light in any direction.
- The gray reference plate is primarily used for optical calibration measurements, such as light source, color, and spectral analysis.
- The grayscale plate has excellent optical stability and never degrades, resistant to acid, alkali, and salt corrosion.

### Specifications

| Product name         | Diffuse Reflection Grayscale Step Plate                                             |
|----------------------|-------------------------------------------------------------------------------------|
| Specification model  | PHL1-S100                                                                           |
| Diameter             | 100mm (30mm-3000mm customizable)                                                    |
| Material             | Spray PTFE                                                                          |
| Spectral range       | 250nm-2500nm                                                                        |
| Reflectivity         | Can be customized according to different requirements for reflectivity              |
| Data and Information | Provide a copy of the reflectance data report from the Chinese Academy of Metrology |
| Shell material\color | Aluminum alloy\black, blue                                                          |

## Spectrum Accessories

### Standard Aluminum Mirror




#### Product Advantages

- The surface of the aluminum mirror is mirror reflection, and the fused quartz coating on the mirror surface prevents oxidation.
- The reflectivity is approximately 85% to 90% at wavelengths ranging from 250nm to 800nm, and 85% to 98% at wavelengths ranging from 800 to 2500nm
- The measurement value is traced back to PRD, and its reflectance ratio is known. A reflectance ratio data report from the National Institute of Metrology, China is provided.

#### Product Description

The standard mirror is a reference material for reflectance measurement. The STA-AL standard mirror from Goptica is made by thermally evaporating high-reflective Al material, capable of providing high reflectance in the wide spectral range of ultraviolet-visible-near infrared (200-2500 nm) for testing reference of mirror reflection samples.

#### Reflectance Spectrum



#### Specifications

| Category     | Parameter                                                                                                         |
|--------------|-------------------------------------------------------------------------------------------------------------------|
| Bandwidth    | 200-2500nm                                                                                                        |
| Material     | Mirror with fused quartz coating                                                                                  |
| Dimensions   | Diameter 38mm, height 11.5mm,<br>Working area aluminum mirror with a 30mm diameter, reflective surface height 8mm |
| Weight       | 32g                                                                                                               |
| Reflectivity | 250-800 nm, ~80%-90%;<br>800-2500 nm, ~85%-98%                                                                    |

## Spectrum Accessories

### Deuterium Halide Tungsten Light Source



#### Product Description


The DHG-L25 composite light source integrates a continuous deuterium lamp and a tungsten-halogen lamp with a broadband spectrum in one channel.

The deuterium lamp emits light primarily used for its ultraviolet wavelength range of 190-400 nm continuous spectral band. The deuterium lamp has three characteristic spectral lines at 486.0 nm, 583.0 nm, and 656.1 nm, which are commonly used as light sources for high-precision analytical measurement instruments, such as those used in liquid chromatographs.

The wavelength range of the tungsten-halogen lamp is typically between 360 nm and 2500 nm. The lifespan of a tungsten-halogen bulb is related to its operating temperature; the higher the color temperature, the shorter the lifespan.

#### Product Features

- It adopts a combination of deuterium lamps and halogen lamps, which has a wider spectral range and can continuously output stable spectra from 190 to 2500nm. The spectral intensity drift is only <0.1%/hr
- Two sets of optical paths are coupled to an SMA905 interface output using a Gewu optical lens. Halogen lamps can be output separately or simultaneously, and the intensity of the output light can be adjusted
- It adopts air-cooled heat dissipation
- It is powered by alternating current ranging from 85 to 250 volts
- The bulb power of the light source is 25W for the deuterium lamp and 10W for the halogen lamp, with a lifespan of 1500 hours
- The size of the light source is 145mm x 200mm x 110mm
- The light source bulb is a consumable. Please replace it before its lifespan expires
- It is recommended to use Gewu Optics' anti-ultraviolet irradiation quartz optical fiber
- Optional external trigger option



#### Product Applications

- Broad spectral range analysis
- Scientific research
- Reflection/Transmission/Absorption spectral analysis
- Deep ultraviolet spectral analysis

## Halogen Light Source




### Product Features

- Adopts high-quality imported halogen bulbs, with the longest lifespan of up to 10,000 hours
- Light source intensity adjustable, adjustment range 0-100%
- Wide wavelength range, 360~2500nm
- Uses air-cooled heat dissipation
- Powered by 85~250V AC or 5V USB directly, convenient for carrying and operation
- Adopts Goptica's lens coupled to SMA905 interface output
- Light source dimensions are small, with the 100W light source measuring 160mm x 85mm x 60mm, and other light sources measuring 150mm x 75mm x 50mm
- The light bulb of the light source is consumable, please replace it before the service life is reached, please contact our after-sales service.

### Product Application

- Color measurement and analysis
- Spectral Analysis of Food and Drug Products
- Spectral Analysis of Environmental Water Pollution



### Specifications

| Model  | Power Consumption | Voltage | Color temperature | Working life | Output optical power | Stability | Dimensions (mm) |
|--------|-------------------|---------|-------------------|--------------|----------------------|-----------|-----------------|
| HG-L5  | 5W                | 5V      | 2800K             | 10000H       | 600um=3.0mW          | <0.15%/hr | 150x75x50       |
| HG-L18 | 18W               | 12V     | 3000K             | 2000H        | 1500um=6.5mW         | <0.2%/hr  | 150x75x50       |
| HG-L9  | 9W                | 5V      | 2900K             | 6000H        | 1000um=5.5mW         | <0.15%/hr | 150x75x50       |
| HG-100 | 100W              | 12V     | 3300K             | 2000H        | 600um=13.0mW         | <0.2%/hr  | 190x92x101      |

## Spectrum Accessories

Standard Light Source Kit:  
Standard A Light Source+Grayscale filter+Cosine corrector

### Standard A Light Source HG-L9A




#### Product Features


- High stability
- Precise light-emitting position
- Multiple accessories available, facilitating various output requirements such as intensity, wavelength, and collimation
- SMA905 output, optional FC output
- High fiber coupling efficiency

#### Product Applications

- Transmissive and reflective measurement
- High-precision illumination



#### Output Spectrum



#### Product Specifications

| Model  | Power Consumption | Power Supply | Color temperature | Working life | Chromatic Coordinate | Interface             | Bandwidth  | Dimensions (mm) | Weight |
|--------|-------------------|--------------|-------------------|--------------|----------------------|-----------------------|------------|-----------------|--------|
| HG-L9A | 9W                | 12VDC 2A     | 2915K             | 6000H        | X:0.444<br>y:0.406   | SMA905<br>or<br>FC/PC | 360-2500nm | 120x55x40       | 290g   |

## Spectrum Accessories

### Standard Light Source Kit:

Standard A Light Source+Grayscale filter+Cosine corrector

#### Grayscale filter



#### Product Description

We can provide grayscale filters with attenuation rates of 10%, 30%, 50% and 70% to achieve the output purpose of gradient luminous intensity. The light source is equipped with a dedicated filter slot.



#### Cosine corrector



#### Product Features

- Wavelength range: 200-2500nm
- SMA905 interface
- Scattering materials: UV-VIS or VIS-NIR

余弦校正器测量太阳辐射



#### Specifications

| Model     | Wavelength range | Dimensions (outer diameter/inner diameter) | Length (external) | Diffuser (diameter) | Diffuser (thickness) | Field of view | Connectable                         |
|-----------|------------------|--------------------------------------------|-------------------|---------------------|----------------------|---------------|-------------------------------------|
| CC-UV-NIR | 200-2500nm       | 6.35mm/4mm                                 | 17.0mm            | 4900μm              | 500μm                | 180°          | Optical fiber with SMA905 connector |

## Spectrum Accessories

### Deuterium Light Source




#### Product Features

- It adopts high-quality bulbs with a lifespan of up to 1,500 hours
- The wavelength range is 190 to 400nm
- The light source output is stable, and the spectral intensity drift is less than 0.1%/hr
- It adopts air-cooled heat dissipation
- It is powered by alternating current ranging from 85 to 250 volts
- The size of the light source is 145mm x 200mm x 110mm
- The light source bulb is a consumable. Please replace it before its lifespan expires
- It is recommended to use Goptica's anti-ultraviolet irradiation quartz optical fiber
- Optional external trigger option

#### Product Applications

- The field of deep ultraviolet spectroscopy analysis
- Ultraviolet absorbance analysis
- Ultraviolet science research
- Chromaticity analysis
- Fluorescence powder analysis



#### Specifications

| Model  | Power consumption | Color temperature | Working life | Interface    | Stability | Weight | Dimension (mm) |
|--------|-------------------|-------------------|--------------|--------------|-----------|--------|----------------|
| DE-L25 | 25W               | 2800K             | 1500H        | SMA905 or FC | <0.15%/hr | 3.5kg  | 145x200x110    |
| DE-L30 | 30W               | 3000K             | 2000H        | SMA905 or FC | <0.2%/hr  | 3.5kg  | 145x200x110    |

## Spectrum Accessories

### Mercury-argon Light Source


#### Product Features



- It adopts high-quality bulbs with a lifespan of up to 6,000 hours
- The size of the light source is 105mm x 75mm x 40mm
- It is directly powered by 12V DC and is convenient to carry and operate
- The Goptica's lens is coupled to the SMA905 interface for output
- The light source bulb is a consumable. Please replace it before its lifespan expires
- It is recommended to use Goptica's anti-ultraviolet irradiation quartz optical fiber

#### Product Applications

- Wavelength calibration light source



#### Specifications

| Model | Voltage | Color temperature | Working life | Interface    | Wavelength range | Weight | Dimension (mm) |
|-------|---------|-------------------|--------------|--------------|------------------|--------|----------------|
| Hg-1  | 12VDC   | 2800K             | 6000H        | SMA905 or FC | 253-1704nm       | 40g    | 105x75x40      |

### Laser



The ThunderBird series of laser sources, excitation sources for fluorescence or Raman research, is characterized by small size and high stability.

- 405nm, line width <5nm, maximum power 100mW
- 488, 532, line width <0.1nm, maximum power 100mW
- 633, 785, 830, 1064nm, line width <0.1nm, maximum power 500mW
- Volume: 45mm x 47mm x 95mm


[www.goptica.com](http://www.goptica.com)

## Spectrum Accessories

### Laser-driven plasma wide-spectrum light source

#### Product Description

Laser-driven plasma white light sources are a new generation of wide-spectrum light sources, mainly applied in many fields such as photolithography, measurement, and spectral detection. This light source uses laser to excite xenon gas to generate plasma light emission. Compared with traditional wide-spectrum light sources, its stability, radiation brightness and lifespan have been improved by several times.



#### Product Features

- High brightness
- Spatial light or optical fiber coupling
- All-in-one computer or discrete computer
- Ultra-wide spectrum
- Long life

#### Product Applications

- Lithography
- Semiconductor measurement
- Spectral detection
- Life sciences and bioimaging

## Spectrum Accessories

### Specifications

| Model                                                             | V199                     | V199FC                   | V29A                     | V29B                     | V277                    | V277H                   |
|-------------------------------------------------------------------|--------------------------|--------------------------|--------------------------|--------------------------|-------------------------|-------------------------|
| Wavelength range                                                  | 170-2100nm               | 170-2100nm               | 170-2100nm               | 170-2100nm               | 170-2100nm              | 170-2100nm              |
| Spectral radiation luminance<br>mW/(mm <sup>2</sup> .sr.nm)@400nm | 20                       | 0.1                      | 10                       | 10                       | 60                      | 100                     |
| Total power                                                       | 0.8W                     | 0.15W                    | 0.5W                     | 0.5W                     | 4W                      | 8W                      |
| Numerical aperture                                                | 0.47<br>(Standard value) | 0.22<br>(Standard value) | 0.56<br>(Standard value) | 0.56<br>(Standard value) | 0.5<br>(Standard value) | 0.5<br>(Standard value) |
| Power stability<br>(RMS)                                          | <0.3%(30S)               | <0.2%(30S)               | <0.2%(30S)               | <0.2%(30S)               | <0.2%(30S)              | <0.2%(30S)              |
| Output method                                                     | Free-space output        | SMA optical fiber output | Free-space output        | Free-space output        | Free-space output       | Free-space output       |
| Service life                                                      | 10000 h                  | 10000 h                  | 10000 h                  | 10000 h                  | 10000 h                 | 10000 h                 |
| Spatial stability of the center of mass<br>RMS                    | <5μm                     | <5μm                     | <5μm                     | <5μm                     | <5μm                    | <5μm                    |
| Luminescent point size(μm)                                        | 100x200                  | 100x200                  | 80x200                   | 80x200                   | 140x450                 | 140x450                 |
| Dimensions of the lamp chamber<br>(mm)                            | 82x86x76                 | 82x86x76                 | 109x120x43               | 109x120x43               | 125x207x95              | 125x207x95              |
| Light source driven ruler(mm)                                     | 113x111x209              | 113x111x209              | none                     | 157x251x133              | 157x251x133             | 157x251x133             |

## Spectrum Accessories

### LED Light Source



#### Product Applications

- Wavelength reference
- Fluorescence measurement
- Lighting

#### Product Features

- All of Goptica's LED light sources are made of imported high-quality lamp beads, featuring stable performance and a lifespan of up to 100,000 hours
- The wavelength is determined by the lamp beads, with options ranging from 245 to 2000nm
- Power supply: 220VAC or 12VDC
- Output interface: SMA905, FC is optional
- TTL trigger is optional

#### LED light source model

| Model      | Electric power | Output optical power | Central wavelength | Half-peak width | Output interface |
|------------|----------------|----------------------|--------------------|-----------------|------------------|
| SQLED-254  | 0.8W           | 0.8mW                | 254nm              | 12nm            | SMA905 or FC     |
| SQLED-275  | 0.8W           | 0.8mW                | 270nm              | 12nm            | SMA905 or FC     |
| SQLED-370  | 3W             | 2~5mW                | 370nm              | 10nm            | SMA905 or FC     |
| SQLED-395  | 3W             | 2~5mW                | 395nm              | 10nm            | SMA905 or FC     |
| SQLED-405  | 3W             | 2~5mW                | 405nm              | 10nm            | SMA905 or FC     |
| SQLED-525  | 3W             | 2~5mW                | 527nm              | 30nm            | SMA905 or FC     |
| SQLED-880  | 2W             | 2~5mW                | 880nm              | 50nm            | SMA905 or FC     |
| SQLED-980  | 2W             | 2~5mW                | 980nm              | 40nm            | SMA905 or FC     |
| SQLED-1550 | 1W             | 1.5mW                | 1550nm             | 120nm           | SMA905 or FC     |
| SQLED-WHT  | 3W             | 2~5mW                | White LED          |                 | SMA905 or FC     |


contact:



+86 150 0085 3620  
sales@goptica.com

## Spectrum Accessories

### Optical Holders AO Series




#### AO-R single-channel TA reflection bracket

- Install a Y-type optical fiber with a TA connector (the connector can be selected by yourself)
- Solid surface reflection measurement
- Excited fluorescence measurement
- Install a cylinder with a diameter of 6.35mm
- The diameter of the fixed hole is 6.35mm



#### AO-RF reflective fluorescence stent

- Install multiple optical fibers (connectors are optional)
- Two excitation light sources can be connected simultaneously
- Fluorescence measurement of solids or powders
- Solid or powder reflection measurement
- Fluorescent filters can be inserted



#### AO-B cuvette holder

- Suitable for 1cm optical path cuvettes
- Measurement of liquid transmittance, absorbance and reflectance
- Enhanced aluminum mirrors can be added
- Applicable wavelengths: 200 to 2500nm
- Two or four SMA905 interface optical fibers are available as options



#### AO-BF cuvette fluorescent scaffold

- Suitable for 1cm optical path cuvettes
- Measurement of liquid transmittance and excitation fluorescence
- Fluorescent filters can be inserted
- Applicable wavelengths: 200 to 2500nm
- Two 180° or 90° SMA905 optical fibers can be optionally connected

## Spectrum Accessories

### Optical Holders AO Series



#### AO-T transmission spectroscopy scaffold

- The length adjustment range of the base plate slide table is 200 to 700mm
- Interface input/output SMA905
- The applicable wavelength range is 200 to 2500nm
- The height of the collimating mirror is adjustable
- A sample rack is optional



AO-TR透反射测试台

#### AO-TR (TRS) reflection and transmission stand

- Solid/liquid reflection measurement
- Solid/liquid transmission/absorption measurement
- Fluorescence measurement
- Scattering measurement
- Diffuse reflection measurement
- Used in conjunction with an optical trap
- 30cm adjustable slide table



AO-TRS透反射积分球测试台

### AO-Dmax optical fiber attenuator

Optical fiber attenuators are used to attenuate the input optical power, avoiding distortion of optical receiving devices caused by extremely strong input optical power.

As a passive optical device, the optical fiber attenuator is used to debug the attenuation of optical signals and adjust the optical power to the required level.




- Physical dimensions: 38\*28\*25mm
- Attenuation wavelength: 200-2500 nm
- Transmittance mode: 0-100%
- Attenuation method: Physical occlusion
- Adjustment method: Knob
- Adjustment method: Fiber interface: SMA905

## Spectrum Accessories

### Fiber Optic Collimator

#### Product Description



Fiber collimating mirrors are used in various optical measurement systems. They are often used to collimate the scattered light output by optical fibers into parallel light, and can also couple parallel light to multimode optical fibers.

The fiber collimating mirrors produced by Goptica's can be used for fiber collimation as well as for coupled focusing. The former can be connected to either single-mode or multi-mode optical fibers. As a coupled focusing, it is only applicable to multimode fibers (numerical aperture 0.22-0.37NA, core diameter  $\geq 100 \mu m$ ).

#### Specifications

| Product model            | ColliM-5                                                                            | ColliM-10                                                                            | ColliM-25                                                                             |
|--------------------------|-------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| Appearance               |  |  |  |
| Lens material            | Ultraviolet fused quartz                                                            | Ultraviolet fused quartz                                                             | Ultraviolet fused quartz                                                              |
| Lens diameter            | 5mm                                                                                 | 10mm                                                                                 | 25.4mm                                                                                |
| Focal length of the lens | 10mm                                                                                | 15mm                                                                                 | 100mm                                                                                 |
| Wavelength range         | 185~2500nm                                                                          | 185~2500nm                                                                           | 185~2500nm                                                                            |
| Interface                | SMA905/FC/TA                                                                        | SMA905                                                                               | SMA905                                                                                |
| Compatible aperture      | 0.22~0.37NA                                                                         | 0.22~0.37NA                                                                          | 0.15~0.22NA                                                                           |
| Shell material           | Aluminum alloy                                                                      | Aluminum alloy                                                                       | Aluminum alloy                                                                        |
| Working temperature      | No frost formation within the temperature range of -20 to 150°C                     | No frost formation within the temperature range of -20 to 150°C                      | No frost formation within the temperature range of -20 to 150°C                       |
| Fixed method             | Thread 3/8-24                                                                       | The outer diameter is 16mm                                                           | Thread M4 or outer diameter 30mm                                                      |